

Electrophilic Aromatic Iodination of Vanillin


Reference: Handout; Chemistry lessons: electrophilic aromatic substitution, recrystallization, melting point;
Green lessons: safer and easier to handle reagents and solvents, more selective reagents

Purpose: To carry out the iodination of vanillin using electrophilic aromatic substitution reactions

Table of Reagents:

Reagents	MW	Amount	BP (°C)	MP (°C)	Density
Vanillin	152.5 g/mol	1 g	285 °C	82 °C	1.056 g/cm ³
Ethanol (EtOH)	46.07 g/mol	20 mL	78.37 °C	-114.1 °C	0.789 g/cm ³
Sodium Hypochlorite	74.44 g/mol	11 mL	101 °C	18 °C	1.11 g/cm ³
Sodium Thiosulfate	158.11 g/mol	10 mL	100 °C	48.3 °C	1.67 g/cm ³
Hydrochloric Acid	36.46 g/mol	~6 mL	48 °C	-26 °C	1.2 g/cm ³
2-Propanol	60.1 g/mol	Varies	82.5 °C	-89 °C	0.786 g/cm ³
Water	18.02 g/mol	Varies	100 °C	0 °C	1.0 g/cm ³
Sodium Iodide	149.89 g/mol	1.17 g	1304 °C	661 °C	3.67 g/cm ³

Balanced Chemical Equation:



Iodination of vanillin

Safety:

- Sodium Hypochlorite → bleach; may damage clothing
- Avoid open flames → methanol/ethanol flammable

Experimental Procedures	Data & Observations
1. Obtain 100 mL round-bottom flask; Add 1 g vanillin with 20 mL ethanol → dissolve; After add 1.17 g sodium iodide, then cool to 0 °C with ice bath	<u>Vanillin Used:</u> 1.0237 g

<p>2. Obtain separatory funnel, add 11 mL of sodium hypochlorite. Add drops to reaction mixture over 10 min period. ***Color will turn from pale-yellow to red brown ***Allow to warm to RT when red brown***</p>	<p><u>Observations:</u> After first 5 drops turned red-brown</p>
<p>3. Add 10 mL sodium Thiosulfate; Add HCl until acidic (use pH paper) ***Add ~6 mL of HCL; precipitate should form***</p>	
<p>4. Remove ethanol from suspension on a rotary evaporator; use gentle heating for <10 min *** Water will still be present***</p>	<p><u>Mass of Product:</u> 0.0223 g <u>Melting Point:</u> 228.6 °C</p>
<p>5. Ice bath the flask for 10 min and collect precipitate by vacuum filtration. →Wash with cold water and ethanol →Continue to dry with vacuum; record mass</p>	<p><u>Crude Product Mass:</u> 0.729 g</p>
<p>6. Recrystallize crude product from ethyl acetate or 2-propanol. Collect crystalline product through vacuum filtration. Weigh product and determine melting point.</p>	<p><u>Pure Product Melting Point:</u> 183 – 185 °C <u>Crude Product Melting Point:</u> 178 °C</p>
<p>7. Recrystallization 2-propanol/Water →Place crude product in 100 mL Erlenmeyer Flask and heat → Add enough 2-propanol to dissolve →Add hot water until mixture cloudy →Add 2-propanol to make clear →Allow to cool and place in ice bath to recrystallize</p>	<p><u>Color:</u> White</p>

Post-lab Questions:

1. a) We collected 0.729 g of the crude product as well as the color of our product was Eggshell-White.

$$\% \text{ Yield} = \frac{\text{Product Obtained}}{\text{Theoretical Yield}} \times 100$$

Obtained: 0.729 g Color: Eggshell-White

- First, we must obtain the theoretical value of our crude product in grams:

$$1.0237 \text{ g Vanillin} \times \frac{1 \text{ mol Vanillin}}{152.15 \text{ g Vanillin}} \times \frac{1 \text{ mol Product}}{1 \text{ mol Vanillin}} \times \frac{278.04 \text{ g Product}}{1 \text{ mol Product}} = 1.871 \text{ g Product}$$

- Next, we must use equation above to get our % Yield:

$$\% \text{ Yield} = \frac{0.729 \text{ g}}{1.871 \text{ g}} \times 100 = 38.96\% \text{ yield}$$

- We received a 38.96% yield for our product
- Our % yield is rather low, which might have resulted in our reagents not mixing properly or has not fully reacted yet


b) Our crude product melted at 178 °C, which concludes that it contains impurities as it has experienced melting point depression. The true melting point range of the pure product was recorded to be between 183 – 185 °C, whereas our product melted below the range. On the same note, this proves that our product's crystalline structure was disrupted by impurities, which lowered our thought to be pure product's melting point (some of the reagents might have been still in the product).

2. On the ^1H NMR we expect to see 5 different resonances of the product, while on the ^{13}C NMR, we would expect to see 8 resonances from the product.

3. $2\text{H}^+ (\text{aq}) + \text{ClO}^- (\text{aq}) + 2 \text{I}^- (\text{aq}) \rightarrow \text{Cl}^- (\text{aq}) + \text{I}_2 (\text{aq}) + \text{H}_2\text{O} (\text{l})$


a) The oxidizing agent is the ClO^-

b) The bleach is environmentally more friendly as it is greener in comparison to using a more powerful oxidizing agent. In addition, generating I_2 produced a mono-iodinated product in comparison to a poly-iodinated product.



a) The practical advantage of reducing excess iodine to iodide is that it is effective and selective in creating the desired regiochemistry. In addition, iodine prefers to stay in a gaseous form as it is a halogen. However, by reducing the excess iodine to iodide we would turn it into an aqueous state, which will allow it to react better with the other reagents. We would overall receive a better yield as the reagents would be able to react better.


5. The HCl reacts with the product in reducing its solubility, which results in a precipitate. Essentially, a H^+ is taken by the polar product in making less soluble



6. $Atom\ Economy = \frac{MW\ of\ products}{MW\ of\ reactants} \times 100$

$$Atom\ Economy = \frac{278.04}{(152.15 + 74.44 + 149.89)} \times 100 = 73.9\%$$

7. The iodination occurs preferentially ortho to the hydroxyl group rather than ortho to the aldehyde ("formyl") group primarily because of resonance stabilization. In this case, para would not be possible in a sense that the hydroxyl and aldehyde groups would be inhibiting that position.



In the instance above, the aldehyde is turned into a carbocation making it a meta director which is a deactivating group.

In this instance, a hydrogen ion is lost in the hydroxyl group resulting in the support for the ortho position. Furthermore, this results in the O^- being a strong activating group, making it a O/P director.

